



Networked Robotics Corp.
825 Chicago Ave, Suite F
Evanston, IL 60202, USA
Toll free: (877) FRZ-TEMP
(877) 379-8367

info@networkedrobotics.com
support@networkedrobotics.com

NetworkedRobotics.com

Holding Magnets for Facilitating Installations

The holding magnets in this package are included to demonstrate an option for wire management both outside and inside many “monitored devices”. They are especially useful in situations where the probe needs to be moved for recalibration or for other reasons. They can be used for standard CATx wire or for that of our digital and analog probes.

The magnets can be used across a wide range of temperature (from -80C to 500C) with only slight loss of magnetization at the extreme temperatures in that range.

In cryofreezers most interiors will be aluminum and therefore not magnetic, but this product can be used at liquid nitrogen temperatures around -196C - however an irreversible slight loss of magnetization of about 10% will occur at these temperatures.

To obtain more of these magnets or for more information

Contact Networked Robotics at sales@NetworkedRobotics.com.